

ZOOM STEREO MICROSCOPE



AB-37A Zoom Stereo Microscope:

- Binocular / Trinocular Viewing Head IPD 55-75mm
- Zoom 1x-4x
- Eyepiece Anti Fungal SWF 10x
- WORKING DISTANCE --- 82mm regardless of magnification.
- TRANSMITTED LIGHT --- 12V10W T-H Lamp / 220V
- INCIDENT LIGHT --- 12V10W T-H Lamp with Reflector / 220V
- Packed In Thermocole Box

Note: Pictures, Specifications & Designs are subject to change without prior notice

ABRON EXPORTS

E-mail No:- info@abronexports.com website <http://abronexports.com>