

PART NO.	Luxeo 4Z	List Price
-----------------	-----------------	-------------------

FAD4144000	ZOOM STEREO MICROSCOPE (BIOLOGY & INSPECTION WORK)	
-------------------	---	--

Stand: Dual light stable stand with Ball-bearing focusing controls

Viewing Bodies: Binocular tube, 45°, IPD 54-74mm

Eyepieces: Wide Field eyepiece WF 10x/20mm focusable paired with eye guards

Zoom ratio: 4.4 : 1

Magnification: 8x-35x

Illumination: Built in LED ring light with +/-30 degree oblique incident illumination settings, Rheostat on ring light provides for intensity control, LED base light for transmitted illumination.

Electrical: Universal power supply (110V-240V, 50/60 Hz) UL, CE approved.

Packed in Styrofoam box, with operation manual, allen wrench, dust cover, cleaning cloth, power cord and spare fuses



**Price
USD 757-00**

**USD
1625-00**



**Digital Trinocular with IUV3100
3 MP camera and software**

**Digital Trinocular with IUV5100
5 MP camera and software**

**USD
2525-00**