



AB-34Bd STEREO ZOOM TRINOCULAR MICROSCOPE WITH DIGITAL CAMERA EYE

- Trinocular Head Angle of viewing: 45°
- Zoom ratio: 1:6.4
- Eye Piece WF10X/20mm
- Objective: 0.7x-4.5x
- Magnification 7x-45x
- Interpupillary distance adjustment range: 55-75mm
- Working distance: 100mm
- Reflected transmitted LED illumination

Model - with Digital Eye Camera (remove eye pice and put it)
 or use in Trino Tube complete with camera + CD Software USB Cable
 can be used on Windows and Mac Book Apple (in VLC Player)
 above with

- AB-34B-D1 Stereo with Camera AB-14a **1.5MP** + software CD for laptop
- AB-34B-D3 Stereo with Camera AB-14a **3MP** + software CD for laptop
- AB-34B-D5 Stereo with Camera AB-14a **5MP** + software CD for laptop
- AB-34B-D6 Stereo with Camera AB-14a **5MP** + software CD for laptop

HIGH RESOLUTION CAMERA