




Centrifuge High Speed 400 RPM AC-330H

Mini Centrifuge - Max. Speed 10000 rpm



Model: Spin-Top

Model: Spin-Top

Micro processor controlled, Micro-Centrifuge machine (Mini Centrifuge) speed 1000 to 10000

Technical Specifications

rpm in step of 1000 rpm. Fitted with 0-15 minute electronic digital timer. Supply completed with 2 standard rotors means

- (i) 6 micro centrifuge of 1.5ml capacity (T-01)
- (ii) 12 micro centrifuge tubes of 0.5ml capacity (T-02)

To work on 220 volts 50Hz AC.

Cat. No.	Description	Packing
T-1	Spin-Top	1

Spare Rotor at extra cost :

Code	Capacity	Packing
T-01	6 x 1.5 ml. tubes	1
T-02	12 x 0.5 ml. tubes	1

Technical Data

Size	Weight	Size	Weight
R.P.M.	RCF (x g)	R.P.M.	RCF (x g)
1000	56	6000	2011
2000	223	7000	2737
3000	503	8000	3575
4000	894	9000	4524
5000	1396	10000	5585

Special Features

- Quiet & reliable
- Small foot print
- Built in power supply
- Digital display of time & RPM
- High quality brushless motor