

VORTEX MIXER



**TEST TUBE
SHAKER
CYCLOMIXER**

H-371 VORTEX MIXER (TEST TUBE VIBRATOR)

Vibration Type
 Variable Speed Vortex Shaker
 Neoprene cup
 in parts vigorous agitation to Test Tube/ Small Flask
 Vibration starts when Test Tube is pressed against cup
 with ON/OFF switch and indicator
 The Cup is having space for Test Tube to prevent walking
 Complete with cord and Plug.
 220V, 50Hz AC operated



Website :

<http://insif.com> abronexports.com

Email sales@insif.com exports@insif.com Ph +91-171-2661618, 4008118 Fax: 2661559 MB 94160-20018

106 MAHESH NAGAR

AMBALA CANTT-133001 INDIA

We Proud to be AN ISO 9001:2000 Company