



WATER BATH

DIGITAL

DOUBLE WALL AC-376



water bath Digital Double walled with Hut Top Lid



with Ring Type Lid 6 holes



with Ring Type Lid 12 holes

SPECIFICATIONS

CONSTRUCTION

Outer Body Double walled Powder Coated

Inner chamber made of Stainless Steel Sheet

Heating By Immersion type heater. 220v 50Hz AC with 2 meter 3 core lead

TEMPERATURE- 5°C above ambient to 100°C X Accuracy 0.5°C

Control- Digital PID Least Count 0.1°C (1/10°C)

LID- Lid with beaker Holes

Model	Sizes	suited for	capacity
AC-376g	270x250x1750mm	6 holes 75mm dia	8 Liters
AC-376h	300x250x175mm	9 holes 75mm dia	12 Liters
AC-376-I	405x355x200mm	12 holes 75mm dia	15 Liters

Optional

AC-376P LID HUT Type **Provided**

AC-376T Test Tube stand 25mm x12 tube x 2 no for water bath

AC-376X Stirrer for water fitting with speed regulator for better accuracy

AC-376Z water circulator motor submircible with speed regulator

Note--Available in HUT Type Lid or Hiles Type Lid any one same price

ABRON EXPORTS

SCIENTIFIC INSTRUMENTS COMPANY

We Proud to be AN ISO 9001:2008 Company

