

U V CABINET



SPECIFICATIONS

Many Specimens of Biology and Geological interest will Fluorescence characteristically when illuminated with U V light

The viewing box allows observation of Bacterial Colonies , Bio-Chemicals , Rocks, etc under U V Light whilst the eyes are shielded from direct radiations.

U.V. lamp with set of 2 of 1 each

- 1 Long wave (366)
- 2 Short wave (254) radiation
- 3 Day light tube

complete wire power supply and wires

U.V viewing cabinet is made of wooden boards with sunmica sheet with glass windows on top and window on side for sample

It also have two openings in front with cloth for inserting hands for sample

Power required as 220v 50hz AC complete with wire and shoe

P-287-05

P-287-01 SUNMICA BOX

P-287-03 WOODEN BOX DELUX MODEL WITH ONE SIDE GLASS WINDOW

P-287-05 PLY BOX WITH ONLY ONE TUBE 25 X 25 X 35 CM



ABRON EXPORTS
A SCIENTIFIC INSTRUMENTS FACTORY



Website : www.insif.com
E-mail: exports@insif.com

AMBALA CANTT-133001 w.e.f. 1-3-2007

We Proud to be AN ISO 9001:2000 Company

13