



B.O.D. ENVIRONMENTAL CHAMBER DIGITAL PID

APPLICATION- **Cat No AC-327-EC**

Water Pollution, Sewerage, Agriculture,
Food Research Lab. Oceanography.

CONSTRUCTION-

Double walled Insulated body
+ having Compresure.

BODY-

Outer chamber Made of Mild Steel Steel,
Powder coated

Inner Chamber Stainless Steel
or Aluminium

Trays generally two perforated trays

Special Door Perxpex
with magnetic Latch for view

Chamber Light CFL Light with auto
off when door is closed

TEMPERATURE Digital PID

-5^o C to 50^o C with cooling arrangement.

ACCURACY- ± 0.5^o C with Digital.

Fan Motor One for cooling.
2nd for Heating.

WORKING- On 220V AC.



Catalogue No	SIZES-	Width	Depth	High	Cap	Cu Ft
AC-327-EC-AA		455	x410	x 610 mm	112 Liter	4 cu.ft.
AC-327-EC-A		505	x 415	x 830 mm.	171 Liter	6.1 cu.ft.
AC-327-EC-B		570	x 550	x 865 mm.	280 Liter	10 cu.ft.
AC-327-EC-C		570	x 550	x 900 mm.	336 Liter	12 cu.ft.
AC-327-EC-D		700	x 650	x 900mm.	420 Liter	15 cu.ft.

OPTIONAL: AC-327-T Timer 0-24 hours 1 day and Night. dial type with notch
AC-327-S Timer 0-24 hours 7 day digital with battery Backup

ABRON EXPORTS
SCIENTIFIC INSTRUMENTS COMPANY

We Proud to be AN ISO 9001:2008 Company

