

Centrifuge High Speed (Brushless Motors)

(Microprocessor Based) 20000 r.p.m.



Model: T-12BL

- High speed Micro Centrifuge,
- Heavy duty C.R.C. construction
- Powder coated finish
- Fitted with brushless automatic balanced motor
- Digital Timer 0-99 minutes
- Digital Speed Indicator
- Digital R.C.F. Indicator
- Programmable Speed regulator
- Electronic Safety lid lock
- Imbalance detector.

Special Features

- Table Top compact
- Brushless self Balancing Motor without carbons
- Ultra High Speed
- Microprocessor based
- Low noise



websites :-
www.abronexports.com

Abron Exports