

**Centrifuge 16000 RPM with timer 0-10m x1m
(Microprocessor Based) also
Model AC-330d2 with timer**

**Model:
AC-330d2**



- Centrifuge Machine, Max. speed 16000 r.p.m
- Compact body which is made of plastic injection moulded.
- Inner safety bowl made of stainless steel
- Fitted with well balance universal motor having long life with imported carbons
- Digital display shows Set time/ Run time /R.P.M.
- Microprocessor based Digital timer of range 0-10 min with LCD display.
- Preset speed controls thru microprocessor from 2000 to 16000 r.p.m.x1000rpm

Supply completed with following accessories:

- Micro rotor to accommodate 8x0.5ml or 1.5ml micro centrifuge

Centrifuge (Brushless) Model: AC-330d2

All specification are same as above but in this model the motor is Brushless (without carbon) **Note:** Speed can also increased upto16000 r.p.m.

Dimension of Machine:	L=290mm, W-270mm, H-240mm
Weight of Machine:	9.5kg. approx.
Dimension of box:	L-355mm x W355 x H-340mm
Weight of Box:	11kg. approx.