

# MINI CENTRIFUGE Abron Cat AB-330HD 10,000RPM Digi



**NEW**

Ideal For Quick Spin-downs  
 Easy To use  
 New Clip-On Designed Rotor  
 Bright LED Display  
 Noise Free Operation  
 Sturdy Compact Design  
 Maintenance Free Motor  
 Large Rotor (helps achieving better Centrifugal force)

Speed : Variable 3000-10000 RPM  
 RCF : 1000g -7500g  
 Timer : 1s-999min  
 Rotor : 8x 0.2ml/0.5ml/1.5ml,2ml  
 2x PCR Tube Strip of 8



Cat No.	Description
MINI-10K+	Mini Centrifuge

\* Cat No. Mini-10K+ Price Includes all Rotors



**Abron Exports** [www.abronexports.com](http://www.abronexports.com)