

**GLASS BEED STERLIZER**  
**MINI FURNACE/ CRUCIBLE FURNACE**

**CAT NO. H-338 GL**

**HANDY  
 &  
 SMALL  
 GIVES VERY HIGH  
 TEMPERATURE  
 IN A  
 SHORT TIME  
 ELECTRICITY SAVING**



**GLASS BEED STERLIZER**

**MINI FURNACE/ CRUCIBLE FURNACE**

- \* MAXIMUM WORKING TEMPERATURE 600-700° C
- \* DOUBLE WALLED ALL STAINLESS STEEL BODY & CHAMBER
- \* WITH MINI DIA 8 INCH( 20 CM ) X 8 INCH HIGH 9 20 CM) 14 LT CHAMBER
- \* HEATERS AROUND CHAMBER
- \* WITH LID FOR BETTER & FASTER TEMPERATURE
- \* FITTED WITH ENERGY REGULATOR
- \* WITH MAIN & LOAD INDICATOR
- \* WORKING ON 220V, 50Hz AC
- \* LOAD 300-500 Watt
- \* FAST & SAFE TO USE
- \* VERSATILE & MULTI USE.

**H-338 GL DIGI**  
**DIGITAL CONTROLLER CUM INDICATOR**  
**0-400'C X 1 'C**



**ABRON EXPORTS**

**A SCIENTIFIC INSTRUMENTS FACTORY**

Website :

<http://insif.com> [abronexports.com](http://abronexports.com)

Email sales@insif.com exports@insif.com Ph +91-171-2661618, 4008118 Fax: 2661559 MB 94160-20018

106 MAHESH NAGAR

AMBALA CANTT-133001 INDIA

**We Proud to be AN ISO 9001:2000 Company**