



Specification	AC-346-HB	AC-343-HA
Temperature Control Range	RT+5°C ~ 100°C	-10~ 100°C
Timer		0 ~ 999 min
Temperature Control Accuracy		±0.5°C
Display Accuracy		0.1°C
Heating Time (20°C To Max Temp.)		≤15min

World's Smallest Dry Bath
 Smallest! Fits in the palm of your hand
 Economically priced
 Easy to Use
 Exchangeable blocks for tubes from 0.2 to 50ml
 Digital temperature control improves temp. Uniformity
 Used advanced thermoelectric technique
 Decent LCD Display with Simultaneous time and temp.
 Beep Sound After program



Cat No.	Description
AC-346HA	Mini Dry bath (RT+5°C to 100°C)
AC-346-HB	Mini Dry Bath (-10to100°C)
AC-346-L	Block For Dry Bath 40 x 0.2ml
AC-346-M	Block For Dry Bath 24 x 0.5ml
AC-346-N	Block For Dry Bath 15 x 1.5ml/2ml
AC-346-O	Block For Dry Bath 4 x 15ml
AC-346-P	Block For Dry Bath 2 x 50ml
AC-346-Q	Block For Dry Bath 8 x 15 ml

