



5893 D&B D-U-N-S® No: 36-908-8502

LAMINAR AIR FLOW BENCH AC-356



Laminar Air Flow



Complete Clean Air Solutions for

- * Biotechnology Units
- * Pathology Labs
- * Pharmaceutical Industries
- * Rresearch Laboratories
- * Institutions



- Made of Stainless Steel/Mild Steel/Ply-board
- Stainless steel working platform
- HEPA filter efficiency 99.97% down to 0.3 micron
- Pre filter efficiency 99% down to 5 micron
- Fluorescent light illumination
- UV germicidal lamp
- Megnahaulic Gauge/Static Pressure Manometer
- Cock for gas/air
- Balanced blower with motor
- Protection grill for HEPA filter
- Removable covers for attending filter and blower assembly
- Micro Processor Controlled unit (Optional)

Working Area Size

- | | |
|---------|-----------------------------|
| AC-356a | Plywood 60x60x60cm |
| AC-356b | Plywood 90x60x60cm |
| AC-356c | Plywood 120x60x60cm |
| AC-356d | Plywood 180x60x60cm |
| AC-356g | Stainless steel 60x60x60cm |
| AC-356h | Stainless steel 90x60x60cm |
| AC-356i | Stainless steel 120x60x60cm |
| AC-356j | Stainless steel 180x60x60cm |

ABRON EXPORTS
SCIENTIFIC INSTRUMENTS COMPANY

We Proud to be AN ISO 9001:2008 Company

