

Sturdy and Compact Design

Super Bright LED display

Variable Speed upto 1500 RPM

Plastic & Ceramic Coated surface are both highly chemical resistant

Microprocessors ensures accurate temperature control

Specifications	AC-359-PS	AC-359-PA
Speed Range	100-1500	100-1500
Max Stirring Capacity	3000ml	3000ml
Dimensions of Plate	5 Inch	5 Inch
Working Plate	Plastic	Ceramic
Display	LED	LED
Heating Output	n/a	500W
Temperature Range	n/a	RT-280°C Increment 1°C



AC-359-PA



AC-359-PS



STIRRING BARS



STIRRING BAR MOVER



Catalog No.	Description
AC-359	FLAT SPIN Magnetic Stirrer
AC-359-PA	Magnetic Stirrer, Plastic Top
AC-359-PS	Magnetic Stirrer + Hot Plate, Ceramic Top
AC-359-S	Stirring bars Assorted Pack 10x6(2Pcs.), 15x8 (2Pcs.), 25x8(2Pcs.), 40x8 (1Pcs.)
AC-359-SB	Stirring Bar Mover