



5893


Magnetic Stirrer Magnetic Stirrer with Hot Plate



AC-359A



BEED




AC-359 B

MAGNETIC STIRRER

- **CONSTRUCTION** :- Mild steel Powder Coated
Top Plate is made of stainless steel.
- **CONTROL** :- Speed regulator is provided
for Controlling Speed of stirrer / beed.
- Hot Plate** :- for heating Solutions
with energy Regulator for heat controll

• MODELS:- All with teflon beed 1pc

- AC-359 A With speed controller
- AC-359 B With speed controller with Hot Plate
- AC-359 C Industrial Model Heavy duty
- AC-359 D Industrial Model Heavy duty with Hot Plate
- AC-359 C With speed controller Industrial Model.
- AC-359 E with speed controll hot plate Digital RPM
- AC-359 H Six in one with speed controll hot plate
all six work independently each have one beed



with Digital RPM AC-359-E



AC-359-C



AC-359H

ABRON EXPORTS

SCIENTIFIC INSTRUMENTS COMPANY

We Proud to be AN ISO 9001:2008 Company

