



FLAT SPIN



STIRRING BARS



STIRRING BAR MOVER

Specifications	AC-359PF	AC-359PA	AC-359PS
Speed Range	15-1500	100-1500	100-1500
Max Stirring Capacity	800ml	3000ml	3000ml
Dimensions of Plate	90x90mm	5 Inch	5 Inch
Working Plate	Plastic	Plastic	Ceramic
Display	n/a	LED	LED
Heating Output	n/a	n/a	500W
Temperature Range	n/a	n/a	RT-280°C Increment 1°C