



Specifications	AC-360D
Speed Range	0-2200
Max Stirring Capacity	20L
Display	LED
Max. Torque	40nm
Overload Protection	Yes
Viscosity Max	10000mPas

High Reliability
 Variable Speed upto 2200 RPM
 Sharp LED Display
 Brushless Motor
 Overload Protection
 Keyless Chuck- Easily Interchange Attachments



Price Includes stand & Stainless Steel Cross Shaft

AC-360d+AC360di



Catalog No.	Description
AC-360D1	Overhead Stirrer (Includes Stand & Stainless Steel Cross Shaft)
AC-360D2	Accessory Straight stirrer, shaft length 40cm, stirrer diameter 6cm, 316L Stainless Steel
AC-360D3	Accessory Blade stirrer, shaft length 40cm, stirrer diameter 6.8cm, 316L stainless steel
AC-360D4	Accessory Centrifugal stirrer, shaft length 40cm, stirrer diameter 9cm, 316L stainless steel
AC-360D5	Accessory Crossed stirrer, shaft length 35cm, stirrer diameter 6.5cm, PTFE coated
AC-360D6	Accessory Straight stirrer, shaft length 35cm, stirrer diameter 7cm, PTFE coated
AC-360D7	Accessory Blade stirrer, shaft length 35cm, stirrer diameter 6.8cm, PTFE coated
AC-360D8	Accessory Centrifugal stirrer, shaft length 35cm, stirrer diameter 8.5cm, PTFE coated
AC-360D9	Accessory Combo Pack - OS-20 SS 6, OS-20 PTFE 6.5, OS-20 PTFE 7