


# Abron Lab Homogenizer AB-360h

## Abron Lab Homogenizer AC-360h

### CONSTRUCTION-

- With fractional H.P. Motor
- working on 220v 50 Hz Operated
- Stainless Steel stirring rod
- Teflon Tip Homogenizer
- with borosilicate glass special tube and holder
- Mounted on heavy cast iron stand,
- with adjustable high clamp
- with speed regulator On / Off

· USAGE-For Tissue homogenization



**Model AC-360k Special Hand high speed**

**ABRON EXPORTS**

SCIENTIFIC INSTRUMENTS COMPANY

Chemistry Heating

We Proud to be AN ISO 9001:2008 Company

Photo & Specifications may change for improvement ( Old Catalogure no was AH-344 ) Website : <http://abronexports.com>