

ABRON DRY HEATING BLOCK DIGITAL CONTROLLED



This Block is comprise a outer case with heater and thermostat, fitted with aluminium blocks drilled to accept test tubes of 25mm dia. Temperature control by means of thermostat. Outside painted in epoxy powder coating outer case is insulated and ventilated.

Temperature range 90 deg. C. To 300 deg C. The above chamber has 40 holes to accept 25 mm test tubes.

SALIENT FEATURES :

1. Aluminium Block Drilled toAccept 25mm dia tube 24 no .
2. Temperature Range 90°C to 300°C with accuracy of $\pm 1^\circ\text{C}$.
3. Digital Temperature controller range cum Indicator