



Shaking Experts

Specifications	AC-366A	AC-366B	AC-366C
Speed Range	50-350 RPM	50-300 RPM	50-250 RPM
Orbit	10mm	20mm	30mm
Timer	1min-999min		
Display	Super Bright LED		
Controller	Programmable Microprocessor		



Compact yet Powerful !
 Extremely Stable and Sturdy even operating at 250 RPM.
 Perfect for Low Capacity Application
 Brushless DC motor ensures smooth operation
 Incubator and Cold room Safe
 Speed and Time Displayed Separate on Large LED Panel
 Choice of Various Platforms



Cat No.	Description
AC-366A	Orbital Shaker 10mm
AC-366B	Orbital Shaker 20mm
AC-366C	Orbital Shaker 30mm
AC-366G	Flat platform with non-slip rubber mat
AC-366H	Universal platform with adjustable bars
AC-366M	Platform with clamps for flasks 100-150ml
AC-366N	Platform with clamps for flasks 250-300ml



**MCOS-FP
AC-366G**



**MCOS-C10250/MCOS-C25300
AC-366M**



**MCOS-UPAB
AC-366H**