

Specifications	AC-H-371M	AC-H-371n
Mixing Motion	Orbital	Orbital
Orbit Diameter	4mm	4.8mm
Speed Range	0-2500	3000
Speed Display	Scale	n/a
Operation	Continuous/Touch	Touch



MX-E

Touch Operation
Fixed 3000 RPM Speed
Small Yet POWERFUL
Eccentric Bearing Design



MX-S

Touch & Continuous Operation
Variable Speed upto 2500 RPM
Heavy Base Avoids unnecessary movement



Accessory Pack For All Tubes 0.2ml, 0.5ml, 1.5ml/2ml, 15ml

Platform Pad

Tube Holding Rod



Cat No.	Description
AC-H-371m	Vortex Mixer Fixed Basic
AC-H-371n	Vortex Mixer Adjustable Speed
AC-H-371-t	Accessories Tube Holding Rod
AC-H-371-s	Accessory Pack Any One (0.2ml, 0.5ml, 1.5ml/2ml, 15ml)
AC-H-371-v	Accessory Pack For All Tubes 0.2ml, 0.5ml, 1.5ml/2ml, 15ml
AC-H-371u	Accessory Pack Flat Platform Pad with Universal Plate