

**CIRCULAR WATER BATH**



**FEATURES H-376-CR      cooling H-376-CCL**

For precise control of different samples at constant temperature. ideal for tissue culture, enzyme reactions fermentations tissue section processing etc. Ideally suited for preparing distortionless and wrinkle free tissue specimens.

Its round, double walled, single piece Inner and outer chambers are made of aluminium, Inside black and outer surface finished in powder paint, with wide enough rim for drying the wet slides. Supplied complete with pilot light, cord and plug.

Suitable to work on 220 V, single phase, 50 Hz AC supply.

Chamber Size :- Dia × Depth × Rim 225 × 70 × 40

- These multisize mantles have heating elements covered in Asbestos/Glass braiding.
- Temperature 40-100° C.
- Built in THERMOSTAT Capillary .
- Lagged with special created mineral wool workable on 220V AC

**WATER BATH WITH COOLING UNIT H-376-CCL**

Temperature -10°c to 30° C.

- Built in cooling unit and heating unit both

Digital temp controller and indicator accuract 0.5°c working on 220V AC



**A SCIENTIFIC INSTRUMENTS FACTORY**

Website : <http://insif.com>  
<http://abronexports.com>

Ph +91-171-2661618, 4008118 Fax: 2661559

Email [sales@insif.com](mailto:sales@insif.com) [exports@insif.com](mailto:exports@insif.com)  
[info@abronexports.com](mailto:info@abronexports.com)

AMBALA CANTT-133001

**We Proud to be AN ISO 9001:2000 Company**

