



# Serological Water Bath Digital PID

Cat. NO. AC-378



Model :- Serological water bath Digital PID Hut Top Lid  
we also have Plain Top Lid

Model :- Powder coated chamber

## SPECIFICATIONS

### CONSTRUCTION

Outer Body Double walled Powder Coated  
Inner chamber made of Stainless Steel Sheet

**Heating** By Immersion type heater.  
220v 50Hz AC

### TEMPERATURE-

5°C above ambient to 100°C X 1°C

**Control-** Digital PID Least Count 1°C

**LID-** \*Provided with top Lid



Water Bath Digital with Stirrer

Model	Sizes	suited for	capacity
AC-378g	300x175x175mm	1 Rack	12 Liters
AC-378h	300x250x175mm	2 Rack	14 Liters
AC-378-I	375x300x175mm	4 Rack	18 Liters
AC-378-j	455x300x175mm	6 Rack	25 Liters
AC-378-k	605x300x175mm	8 Rack	32 Liters
with STIRRER	Sizes	Suited for	capacity
AC-378-is	375x300x175mm	4 Rack	18 Liters
AC-378-js	455x300x175mm	6 Rack	25 Liters

**ABRON EXPORTS**  
SCIENTIFIC INSTRUMENTS COMPANY

We Proud to be AN ISO 9001:2008 Company

