


# Balance Digital

210gm 0.1mg Abron AC-255d3r



## **digital balance 30gm/1mg:**

1. Pan: Circle 80 mm (SST 304), Outline Size: 34\*21\*17cm
2. Display: LCD (white back-light) (LxW) 110x30 mm
3. Interface: Rs232C
4. Power: AC (220-240V 50HZ / 110-120V 60HZ)
5. Units: g, oz, ct, lb, GN etc.
6. More functions: counting, last remember, %, full range tare, over-load protection, zero tracking, low-battery indicator etc.

## **Packing:**

- 1Pc / carton
- 42x31x43cm / carton
- G.W.: about 7kg / carton