

**ABRON EXPORTS**  
Formally :- Indian Scientific Instruments Factory

106, Mahesh Nagar,  
Ambala Cantt. - 133 001 (India)

Ph. : 0171- 4008118

• FaxNo. 0171-2661559 • [www.abronexports.com](http://www.abronexports.com) • E-mail : [info@abronexports.com](mailto:info@abronexports.com)

Exporter & Importer of  
Scientific Instruments

TIN No. : 08021007095

Export IEC No 3394003350.

ESTD. 1958 (49 years old)  
Golden Jubilee in 2008

**Abron**  
Instruments

## ABRON BIO FREEZE



## BIO FREEZE

Capacity: 90L, 170L, 290L, 350L, 450L, 600L.

Models: Vertical & Horizontal.

Temperature: -20°C and - 40° C.

Microcontroller based controller.

Inner chamber SS 304, SS trays.

Audio visual alarm. PID Digital controller

**ABRON EXPORTS** India [www.abronexports.com](http://www.abronexports.com)