

# ABRON EXPORTS

Factory :-105 Industrial Area **AMBALA CANTT. 133 001**

NEW Telephone No. +91-171-2699478 2698459 Mobile +91-94160-20018

E-mail No:- [insif@sancharnet.in](mailto:insif@sancharnet.in) [exports@insif.com](mailto:exports@insif.com)

## Autoclaves - M.O.U Portable STERILIZER INSTRUMENT COOKER TYPE

Electric Type  
Model AB-E-SS-10  
Model AB-E-AL-11



The Portable Single Drum Autoclave is A highly efficient apparatus which is designed to meet the requirements of small dispensaries, nursing homes, laboratories, medical missions, clinics, ships operating theatres, pharmaceuticals, etc. It is very simple to operate and economical to use. The container is seamless made of thick Aluminium sheet or stainless steel as per quality required. The lid is made of alluminium casting or stainless steel as per customer requirement. It is fitted with timer, low water cut off pressure switch or any other customer requirements. The safety valve is adjusted to blow off the steam at 17 lbs., p.s.i. maintaining the constant pressure of 15 lbs. The entire autoclave is heavily Aluminium alloy/stainless steel. The container being seamless, there is no leakage and hence it ensures quite a longer service. The autoclave is made for use on electrically and gas.

**These are offered in the following quantities:** 9 lt, 20 lt, 24 lt, 33 lt, 39 lt

**These are available in the following types:**

- Nut Type
- P Type (cooker type)
- Vertical Autoclave



Non Electric Stove Type/  
Kerosene Oil Type  
Model AB-K-SS-12

**Ask for Quote**

These sizes are available in seamless & jointed range. Besides this, we also offer as per our customers' specifications, as per their requirements of quality & range.