



5893

No: 86-908-6302

ABRON AUTOCLAVES Tripple walled Thermostatic

AC-325-D Tripple Walled Abron Vertical Autoclave All Stainles steel with pressure Guage

with safety valve, Pressure Gauge
Automatic Pressure Release Valve
Locking Screws top lid
and water Level Teflon Glass Window
Outer + inner **Body**
Heavy stainless steel sheet
With Stainless steel wire Basket
Three Legs stand



Chamber Dia X Hight inch- (CM)

AC-325D-10	Chamber 10" ht. 18"	250x450mm
AC-325D-12	Chamber 12" ht. 20"	300x500mm
AC-325D-14	Chamber 14" ht. 22"	350x550mm
AC-325D-16	Chamber 16" ht. 24"	400x600mm

Chamber Dia X Hight inch- (CM)

with Digital Programable Controller

AC-325DD-10	Chamber 10" ht. 18"	250x450mm
AC-325DD-12	Chamber 12" ht. 20"	300x500mm
AC-325DD-14	Chamber 14" ht. 22"	350x550mm
AC-325DD-16	Chamber 16" ht. 24"	400x600mm

OPTIONAL Digital Pressure controler

AC-325DDP Digital Pressure controller
in any one above

Lid Parts



Water Level

AC-325-R is Our new Catalogue Model of AH-376R

ABRON EXPORTS
SCIENTIFIC INSTRUMENTS COMPANY

We Proud to be AN ISO 9001:2008 Company

