

**Abron Eddy's Hot / Cold Plate New Model
APH-2445**



For determining Analgesic effect of morphine in case of mice.

Consisting of heating surface with perpex enclosure with Micro-controller based Digital Temperature Indicator controller to set surface temp.

Between 3C to 80°C.

Price : Rs 0 each

[Send Enquiry](#)

[Modify](#) | [Delete](#)

ABRON EXPORTS 105 Industrial Estate AMBALA CANTT 133001

Ph 0171-2699478 MB 94160-20018 Email info@abron.in www.abronexports.com