



AB-05a MEDICAL MICROSCOPE Abron

Triple Revalving nose
Fixed condenser Iris Diaphragm **Focusable**
Sturdy body, Horse Shoe base.
Coarse & Fine focussing knobs.
Stage clips, Reflector, Dust cover, Cloth.
Coaxial -Mechanical Stage,

Eye piece	Objective	Magnification	
Hygenion	10x,	100x ,	150x
10x, 15x	40x	400x,	600x

AB-5B MEDICAL MICROSCOPE 100X

-DO- with 100x oil immertion lens

Eye piece	Objective	Magnification	
Hygenion	10x,	100x ,	150x
10x, 15x	40x	400x,	600x
	100x	1000x,	1500x
	Oil Immersion Lens		

OPTIONAL

AB-57B

LED Light
with 2x1.5v Cell



AB-57C

LED Light
Variable
2 x 1.5v Cell

AB-57D

LED Light
220v Operated

