

**Balance Digital Lab Table / Analytical
Round air shield with power adaptor**

Plastic with Air Shield + power adaptor 220v

<u>Cat No</u>	<u>Type</u>	<u>Capacity</u>	<u>L C</u>
AC-255c2	Balance digital 200gm x 1mg round air shield		
AC-255c3	Balance digital 300gm x 1mg round air shield		
AC-255c5	Balance digital 500gm x 1mg round air shield		



**Balance Digital Lab Table / Analytical
Air shield with power adaptor**

Plastic with Air Shield + power adaptor 220v



<u>Cat No</u>	<u>Type</u>	<u>Capacity</u>	<u>L C</u>
AC-255c2	Balance digital	200gm x 1mg	air shield
AC-255c3	Balance digital	300gm x 1mg	air shield
AC-255c5	Balance digital	500gm x 1mg	air shield

**Balance Digital Analytical Precision 0.1mg
4th Fourth Digit 220gm x 0.1mg AC-255D4**

Rectangular with
Maximum Capacity
Multiple weighing
Response Time
Printer Output
Shape of PAN
Power Supply

Air Shield + power adaptor 220v
220 g x 0.1mg
Units also piece counting
in 0.2 second
RS232C Interface
Concave
Single Phase



<u>Cat No</u>	<u>Type</u>	<u>Capacity</u>	<u>L C</u>
AC-255D4E	Balance digital 220gm x 0.1mg External calibration	RS232	
AC-255D4i	Balance digital 220gm x 0.1mg Internal Calibration	RS232	
AC-255D4W	Balance digital 220gm x 0.1mg (0.0001) Wireless sensor		

**Balance Digital Micro Analytical High Precision
200gm 0.01mg 5th Fifth Digit AC-255D5E**

Rectangular with
Maximum Capacity
Multiple weighing
Response Time
Printer Output
Date time ID Setting
PAN Concave

Air Shield + power adaptor 220v
220 g
Units also piece counting
in 0.2 second
RS232C Interface
Price counting % percentage setting
Power Single Phase 220v AC adaptor

<u>Cat No</u>	<u>Type</u>	<u>Capacity</u>	<u>L C</u>
AC-255D5E	Balance digital Precision 220gm x 0.01mg RS232 External calibration Date time ID Setting Price counting % percentage setting	Battery Backup	

