

# Binocular Research Microscope LED Light +100x Oil Immersion Lens



**AB-11 DIGITAL BINOCULAR RESEARCH MICROSCOPE WITH CAMERA USB ABRON**

**Focussing-** for coarse/fine  
**Mechanical Stage-**145x 125mm  
**Binocular Head-**High Resoulution  
 Inclined Rotatable 360°  
 A built-in base 6V 20W halogen light  
 with intensity regulator.

**EYE PIECES-** Wide field 10x paired 15x Hygenion.

**OBJECTIVE-** 4x, 10x 40x,  
 100x Oil immersion

with  
**AB-12DC CAMERA 5MP USB CABLE**  
 + software CD  
 for laptop Use

