

Micropipette Digital Variable Volume

Two Button Aircushion Micropipette

Abron AC-289V



- Variable with Eject Tips Easy setting by simply turning mechanical
- Click Stop System for Aqueous Liquids Two Button operation
- Easy to calibrate and service. Pipetting of aqueous liquids,
- two button operation Volume range as below
- Available options: 1-channel Autoclavable for lab uses &
- Purity grades of consumables: PCR clean and DNA forensic
- Autoclavable tip cone and lower section
- Accurate * Conform to DIN 12650 Standards

| Cat No | Type | Volume ml | Least count LC |
|-------------|--------------------------------|-------------|----------------|
| AC-289-V10 | Micro Pipette Digital variable | 0.5-10 ul | 0.1 ul |
| AC-289-V20 | Micro Pipette Digital variable | 2-20 ul | 0.1 ul |
| AC-289-V50 | Micro Pipette Digital variable | 5-50 ul | 1 ul |
| AC-289-V1H | Micro Pipette Digital variable | 20-100 ul | 1 ul |
| AC-289-V2H | Micro Pipette Digital variable | 20-200 ul | 1 ul |
| AC-289-V1T | Micro Pipette Digital variable | 200-1000 ul | 10 ul |
| AC-289-V1TA | Micro Pipette Digital variable | 100-1000 ul | 10 ul |
| AC-289-V5T | Micro Pipette Digital variable | 1-5 ml | 100 ul |
| AC-289-V10T | Micro Pipette Digital variable | 2-10ml | 100 ul |