

## Microscope

### Microscope Student Objective 10x 45x

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X	100 X, 150x
15X	45 X	450X , 675X

**Viewing Head** : Monocular 90°  
**Nosepiece** : Triple revolving Upright  
**Objective** : Achromatic 10x, 40x Spring Loaded SL  
**Stage** : 110x110mm with stage clip  
**Focusing Knobs** : Coarse and Fine separate  
**Condenser** : **No just disc Diaphragm**  
**Illumination** : Plano-concave mirror  
**Carrying Case** Styrofoam **Supplied with:** Dust cover

**Cat No**      **Type**      **Optics**

AB-01-A    Microscope student disc diaphragm,  
Eye piece 10x,15x Objective 10x,45x

AB-01-B    **Microscope student Iris diaphragm,**  
Eye piece 10x,15x Objective 10x,45x

AB-01-B2 **Microscope student Iris diaphragm,**  
Eye piece 10x,15x Objective 10x,45x 100x, oil immersion

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X, 45x	100 X, 150x, 450x
15X	100x oil immersion	675X, 1000x, 1500x



**Abron Exports**  
www.abronexports.com

### Microscope Student 360 Rotatable Tube Ab-02-B

**Cat No**      **Type**      **Optics**

AB-02-B    Microscope student disc diaphragm, 360° Tube  
Eye piece 10x,15x Objective 10x, 45x

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X	100 X, 150x
15X	45 X	450X , 675X

**Viewing Head** : Monocular 360° Rotatable  
**Nosepiece** : Triple revolving Upright  
**Objective** : Achromatic 10x, 40x Spring Loaded SL  
**Stage** : 110x110mm with stage clip  
**Focusing Knobs** : Coarse and Fine separate  
**Condenser** : **Iris Diaphragm**  
**Illumination** : Plano-concave mirror  
**Carrying Case** Styrofoam **Supplied with:** Dust cover

**Optional Accessories:**

AB-57-C      LED Lamp, Variable 2x1.5v Cell  
 AB-54-C      **Eye Piece 10X WF**



**Abron Exports**  
www.abronexports.com

## Microscope Digital

### Microscope Student with Digital Camera Eye

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X	100 X, 150x
15X	45 X	450X , 675X

**Viewing Head** : Monocular 90°  
**Nosepiece** : Triple revolving Upright  
**Objective** : Achromatic 10x, 40x Spring Loaded SL  
**Stage** : 110x110mm with stage clip  
**Focusing Knobs** : Coarse and Fine separate  
**Condenser** : **No just disc Diaphragm**  
**Illumination** : Plano-concave mirror  
**Carrying Case** Styrofoam **Supplied with:** Dust cover

**Illumination** : **LED light & Plano-concave mirror**

**Cat No**      **with software for laptop**      **Mega Pixel**

AB-01-D1    Microscope student with Digital camera software 1.5 MP  
Eye piece 10x,15x Objective 10x,45x

AB-01-D3    Microscope student with Digital camera software 3.0 MP  
Eye piece 10x,15x Objective 10x,45x

AB-01-D5    Microscope student with Digital camera software 5.0 MP  
Eye piece 10x,15x Objective 10x,45x

