

Micro pipette Digital Variable Super Adjustable Volume Micropipettes

Variable with Eject Tips Easy setting by simply turning
 * Click Stop Digital System * Easy to calibrate and service.
 * Autoclavable tip cone and lower section
 * Accurate * Conform to DIN 12650 Standards

Cat No	Type	Volume ml	Least count	LC
AC-289-S10	Micro Pipette Digital variable	0.5-10 ul	0.1 ul	
AC-289-S20	Micro Pipette Digital variable	2-20 ul	0.1 ul	
AC-289-S50	Micro Pipette Digital variable	5-50 ul	1 ul	
AC-289-S1H	Micro Pipette Digital variable	20-100 ul	1 ul	
AC-289-S2H	Micro Pipette Digital variable	20-200 ul	1 ul	
AC-289-S1T	Micro Pipette Digital variable	200-1000 ul	10 ul	
AC-289-S5T	Micro Pipette Digital variable	1-5 ml	100 ul	
AC-289-S10T	Micro Pipette Digital variable	2-10ml	100 ul	



Micro pipette 8 Channel Digital Variable Adjustable Volume Semi Autoclavable

* Variable with Eject Tips Easy Volume setting by turning
 * Click Stop Digital System* Easy to calibrate and service.
 * Autoclavable tip cone and lower section
 * Accurate and precision guaranteed
 * Conform to DIN 12650

Cat No	Type	Variable	Volume ml	Least count	LC
AC-289-M81	Micro Pipette 8 Chanel Digi vari	0.5-10 ul	0.1 ul		
AC-289-M82	Micro Pipette 8 Chanel Digi vari	2-20 ul	0.1 ul		
AC-289-M85	Micro Pipette 8 Chanel Digi vari	5-50 ul	1 ul		
AC-289-M810	Micro Pipette 8 Chanel Digi vari	10-100 ul	1 ul		
AC-289-M820	Micro Pipette 8 Chanel Digi vari	20-200 ul	1 ul		
AC-289-M830	Micro Pipette 8 Chanel Digi vari	30-300 ul	10 ul		



Micro pipette 4 Channel Digital Variable Adjustable Volume Semi Autoclavable

Cat No	Type	Variable	Volume ml	Least count	LC
AC-289-M41	Micro Pipette 4 Chanel Digi vari	0.5-10 ul	0.1 ul		
AC-289-M42	Micro Pipette 4 Chanel Digi vari	2-20 ul	0.1 ul		
AC-289-M45	Micro Pipette 4 Chanel Digi vari	5-50 ul	1 ul		
AC-289-M410	Micro Pipette 4 Chanel Digi vari	10-100 ul	1 ul		
AC-289-M420	Micro Pipette 4 Chanel Digi vari	20-200 ul	1 ul		
AC-289-M430	Micro Pipette 4 Chanel Digi vari	30-300 ul	10 ul		

Micro pipette 12 Channel Digital Variable Adjustable Volume Semi Autoclavable

Cat No	Type	Variable	Volume ml	Least count	LC
AC-289-M121	Micro Pipette 12 Chanel Digi vari	0.5-10 ul	0.1 ul		
AC-289-M122	Micro Pipette 12 Chanel Digi vari	2-20 ul	0.1 ul		
AC-289-M125	Micro Pipette 12 Chanel Digi vari	5-50 ul	1 ul		
AC-289-M126	Micro Pipette 12 Chanel Digi vari	10-100 ul	1 ul		
AC-289-M127	Micro Pipette 12 Chanel Digi vari	20-200 ul	1 ul		
AC-289-M128	Micro Pipette 12 Chanel Digi vari	30-300 ul	10 ul		



Micro pipette Digital Variable Electronic Adjustable Volume Micropipettes

Variable with Eject Tips Easy setting by push bottom

- Based on Stepper Motor Technology
- user friendly software
- Recharging during operation
- Optimized ergonomics and working comfort
- Built-in tip ejector Pipetting range from 0.2 µl to 10000 µl

Cat No	Type	Variable	Volume ml	Least count	LC
AC-289-E1	Micro Pipette Digital variable	0.5-10 ul	0.1 ul		
AC-289-E5	Micro Pipette Digital variable	5-120 ul	0.5 ul		
AC-289-E20	Micro Pipette Digital variable	20-300 ul	1 ul		
AC-289-E50	Micro Pipette Digital variable	50-1000 ul	5 ul		
AC-289-E1H	Micro Pipette Digital variable	100-5000 ul	10 ul		
AC-289-E5H	Micro Pipette Digital variable	500-10000 ul	50 ul		

