



AB-11 DIGITAL BINOCULAR RESEARCH MICROSCOPE WITH CAMERA USB ABRON

Focussing- for coarse/fine

Mechanical Stage-145x 125mm

Binocular Head-High Resoulation

Inclined Rotatable 360°

A built-in base 6V 20W halogen light with intensity regulator.

EYE PIECES- Wide field 10x paired 15x Hygenion.

OBJECTIVE- 4x, 10x 40x, 100x Oil immersion

with

AB-12DC CAMERA 5MP USB CABLE

+ software CD

for laptop Use

