



OPTIONAL

AB-57B LED Light
with 2x1.5v Cell



AB-57D LED Light
220v Operated
with wire plug



AB-57C LED Light
Variable
2 x 1.5v Cell



AB-04A JUNIOR MEDICAL MICROSCOPE STAGE CLIP+ ADJUSTABLE CONDENSER LENS

- Triple Revalving nose with Abbe condense NA 1.2
- Sturdy body, Horse Shoe base.
- Coarse & Fine focussing knobs.
- Stage clips, Reflector, Dust cover, Cloth.

Huygenion Eye Piece 10x and 15x

Achromatic Objectives 10x and 45x

MAGNIFICATIONS 100X 150X 450X, 675X

AB-04b JUNIOR MEDICAL MICROSCOPE+LED OPTIONAL

AB-04d1 Microscope +Digital 1.5 MP Camera

AB-04d2 Microscope +Digital 3 MP Camera

AB-04d3 Microscope +Digital 5 MP Camera

