


2021-23

**MICRO PIPETTE GLASS** Please see Pipette

**Micro Pipette Graduated Serological Boro Glass also catalogue No AC-2212**

Screen printed Suitable for Micro Biology / lab uses  
Graduated to lip and calibrated to Deliver their total capacity when the last drop is blown out

Cat No	Type	Cap ml	LC
AC-457-G01	Micro Pipette Glass Graduated	1 ml	0.01 ml
AC-457-G02	Micro Pipette Glass Graduated	2 ml	0.02 ml
AC-457-G05	Micro Pipette Glass Graduated	5 ml	0.05 ml
AC-457-G10	Micro Pipette Glass Graduated	10 ml	0.1 ml
AC-457-G20	Micro Pipette Glass Graduated	20 ml	0.5 ml
AC-457-G25	Micro Pipette Glass Graduated	25 ml	0.5 ml
AC-457-G50	Micro Pipette Glass Graduated	50 ml	1 ml
AC-457-G1H	Micro Pipette Glass Graduated	100 ml	1 ml



Abron Exports  
www.abronexports.com

**Micro Pipette Graduated Serological Sterile Disposable polystyrene**

Individually Wrapped, Accurate and Bold graduations, Color-coded ring for easy identification, Non Pyrogenic, Dnase Rnase Free Gamma Irradiation Sterilized

Cat No	Type	graduated coded	Cap ml	LC
AC-532-PS01	Micro Pipette polystyrene sterile		1 ml	0.01 ml
AC-532-PS05	Micro Pipette polystyrene sterile		5 ml	0.05 ml
AC-532-PS10	Micro Pipette polystyrene sterile		10 ml	0.1 ml
AC-532-PS25	Micro Pipette polystyrene sterile		20 ml	0.5 ml
AC-532-PS50	Micro Pipette polystyrene sterile		25 ml	0.5 ml




Abron Exports  
www.abronexports.com

**Micro Pipette Fixed Volume Mini AC-289N**

Fixed Volume \* Fully Autoclavable accuracy 1ul

Cat No	Type	Volume ml accuracy
AC-289-N	Micro Pipette Fixed volume Autoclavable	5ul, 10ul, 20ul, 25ul, 50ul, 100ul, 200ul, 250ul, 500ul, 1000ul



**Micro Pipette Digital Fixed Volume AC-289F**

Fixed Volume with Eject Tips Accurate

\* Autoclavable tip cone and lower section

Cat No	Type	Volume ml	Accuracy 1ul
AC-289-F05	Micro Pipette Fixed	5 ul	0.5 ul
AC-289-F10	Micro Pipette Fixed	10 ul	0.5 ul
AC-289-F20	Micro Pipette Fixed	20 ul	0.5 ul
AC-289-F25	Micro Pipette Fixed	25 ul	0. ul
AC-289-F50	Micro Pipette Fixed	50uL	1 ul
AC-289-F1H	Micro Pipette Fixed	100uL	1 ul
AC-289-F2H	Micro Pipette Fixed	200uL	1 ul
AC-289-F25H	Micro Pipette Fixed	250uL	10 ul
AC-289-F1T	Micro Pipette Fixed	1000uL	10 ul
AC-289-F5T	Micro Pipette Fixed	5000uL	100 ul
AC-289-F10T	Micro Pipette Fixed	10000uL	100 ul



Abron Exports  
www.abronexports.com

**Micro pipette Digital Variable Adjustable Volume Micropipettes**

Variable 0-100/200/500/1000 ul with Eject Tips

- \* Easy Volume setting by simply turning
- \* Click Stop Digital System
- \* Easy to calibrate and service.
- \* Autoclavable tip cone and lower section
- \* Accurate and precision guaranteed \* Conform to DIN 12650

Cat No	Type	Volume ml	Least count	LC
AC-289-A10	Micro Pipette Digital variable	0.5-10 ul	0.1 ul	
AC-289-A20	Micro Pipette Digital variable	2-20 ul	0.1 ul	
AC-289-A50	Micro Pipette Digital variable	5-50 ul	1 ul	
AC-289-A1H	Micro Pipette Digital variable	20-100 ul	1 ul	
AC-289-A2H	Micro Pipette Digital variable	20-200 ul	1 ul	
AC-289-A1T	Micro Pipette Digital variable	200-1000 ul	10 ul	
AC-289-A5T	Micro Pipette Digital variable	1-5 ml	100 ul	
AC-289-A10T	Micro Pipette Digital variable	2-10 ml	100 ul	



Abron Exports  
www.abronexports.com