

STUDENT MICROSCOPE ABRON SM-11-LED

with LED Main 220v Light

- Triple Revalving nose with Abbe condense NA 1.2
- Sturdy body, Horse Shoe base.
- Coarse & Fine focussing knobs.
- Stage clips, Reflector, Dust cover, Cloth.

HUYGENIEN EYE PIECE

10x and 15x

ACHROMATIC OBJECTIVES

10x and 45x

MAGNIFICATIONS 100X 150X
450X, 675X



OPTIONAL

AB-57B LED Light
with 2x1.5v Cell



AB-57C LED Light
Variable
2 x 1.5v Cell

AB-57D LED Light
220v Operated
with wire plug

