

Millipore Filtration set for 47mm paper + Vacuum Pump setup Abron AC-471



Having Borosilicate glass Parts

- A Flask 1000ml B Filter 47mm Dia paper Holder with disc
C Buckner funnel 300ml D Filter paper holder Aluminum clamp

Cat No	Type	Size Shape	Unit
AC-471-MF	Filter Flask 1000ml	Holder Borosilicate glass	setup
AC-471-MK	Millipore Filter Flask	1000ml holder	vacuum Pump set
AC-607-M45	Filter Paper Millipore	45µm Dia 47mm	Pack of 100 Pc
AC-607-M22	Filter Paper Millipore	22µm Dia 47mm	Pack of 100 Pc

Vacuum Pump Single Stage Mono block with Oil Nozel for Air 25 Litre / Minutes AC-288A25

Cat No Item Type Air Capacity
AC-288A25 Vacuum Pump Electric 220v
+Oil mono block 25 Lt/ M
Free Air (CFM) 2.5 Vacuum (mbar) 0.003
Displacement (L/min) 25 Litre/min or 50 Litres /
Minutes
Ultimate (pa) 5Pa 0.05mbar
Motor ¼ HP Voltage (CFM)220v 50Hz AC
Intake fitting 1/4" Flare Oil capacity 280-300ml



Abron Exports
www.abronexports.com