


Labomed Digital Camera iVu3100 3mp (Abron Cat no AB-37g 31)



Dimensions (Units: mm)



Sensor	1/2" 3.0 MP Color CMOS
Active area	4.52mm x 3.40mm
Sensor resolution	2048 x 1536 pixels
Pixel Size	2.2µm x 2.2µm
A/D conversion	12 Bit RGB
Pixel clock	48 MHz
Dynamic range	56 dB
Max. exposure	2 seconds (Automatic Exp, WB, BL)
Frame rate (image size)	30 fps (640 x 480)
Image resolution (Standard Binning)	2048 x 1536 1024 x 768 640 x 480
Digital interface	USB 2.0
Optical connection	C- Mount (0.63x)
Voltage supply	5VDC (via USB 2.0 connector)
Power consumption	approx. 1.9W
Ambient conditions	Temperature: -30°C... +70°C Humidity: 5%... 80%, not condensing
Dimensions (LxWxH)	145mm x 93mm x 123mm
Weight	600g
Capture software	PixelPro
Computer Requirements	PC: Microsoft Windows® XP/ Vista/ 7 OS, 3 Ghz CPU, 2GB RAM (if 1GB, then 64MB graphics card recommended), USB 2.0

Package Contents <ul style="list-style-type: none"> • CapturePro Image Capture Software • Firewire IEEE 1394 Cable (3ft) 	Warranty Information 2 Years Camera
---	---