

# Binocular Research Microscope LED Light +100x Oil Immersion AB-11bb with Battery Backup



## AB-11BB BINOCULAR RESEARCH MICROSCOPE ABRON with Battery Backup

**Focussing-** for coarse/fine

**Mechanical Stage-**145x 125mm coaxial

**Binocular Head-High Resoulution**

Inclined Rotatable 360°

**A built-in base LED with variable intensity with battery**

**FILTER** Blue and Green Filter

**EYE PIECES-** Wide field 10x paired  
15x Hygenion. PAIRED

**OBJECTIVE-** 4x, 10x 40x,  
100x Oil immersion

**MAGNIFICATION -** 40x, 60X  
100x 150X 400x, 600X, 1000X, 1500X

