


Optical Bench Single Rod

AP-839-A Optical Bench set Single Rod 1 meter long with 4 riders 2 pins 2 lens holder complete

Single Rod One meter long square rod 18 mm x 18 mm deeply graduated in cm mm two side cast iron stands with one side levelling screw with four sliders Two of the sliders are provided with lateral motion. Complete with two lens holders and two index pins



AP-839-M Optical Bench Single Aluminium pipe 150cm with lamp house 4 riders 2 pins 2 lens holder 12V DC Bulb complete set

Single Aluminium pipe rectangular 150cm long with scale cm mm with die cast side feet with leveling screws with six riders accessories with stem 1 lens holder and 1 index pin, 2 screens, 1 pin round hole screen, 1 light set halogen bulb 12vdc



AP-839M5 optical bench single Aluminium pipe with scale 150cm with riders

AP-839-P Optical Bench Single Rod Spare parts 4 riders 2 pins 2 lens holder complete set



Spare 2 sets Lens holder Pins



Power supply 12v in steps



AP-839-P Optical Bench Single Rod Spare parts 4 riders 2 pins 2 lens holder complete set graduated in cm mm



AP-839-M Optical Bench Single Rod parts 10 riders 2 pins 2 lens holder screen, slit, cross, lamp 12vdc complete set graduated



Abron Exports
www.abronexports.com