

Binocular Research Microscope LED Light +100x Oil Immersion AB-11E



AB-11 BINOCULAR RESEARCH MICROSCOPE BALL BEARING NOSE PIECE CO-AXIAL-STAGE

Focussing-	for coarse and fine focusing Knobs
Mechanical Stage-	145x 125mm XY Rack & pinion co-axial
Binocular Head-	High Resolution Inclined 45° Rotatable 360° inter-pupilear distance adjustable
Condenser	Abbe condenser adjustable N.A. 1.25
Light Variable in base	6V 20W halogen light 100-220v 50Hz/LED
Solar Reflector	for day light work reflector Plano-concave
EYE PIECES-	Wide field 10x paired anti-fungal Blue coated 16mm
Optional	also 15x Hygenion.Eye Piece
OBJECTIVE-	Achromatic Anti-fungal 4x, 10x 40x, 100x Oil immersion Lens
Magnification	40x, 100x 400x 1000X with 15X Eye piece 60x 150x 600X, 1500X
Accessories (1 Eh)	Dust Cover , cleaning Kit, Power cord(tandem Plug), Blue filter , 8ml Oil lens , English Manual
Optional	Wooden Box , Spare Bulb/LED, Spare Filter

