

A b

Electrode Combined for pH meter

With BNC cable connector and semipermeable membrane tip allows electrons to pass Dip in solution bottles

Cat NO	Type dia	Size Shape	Unit
AC-308-P	Electrode combined for pH meter		
AC-308-PC	Electrode Coleman for pH meter		



Electrode Combined stand for pH meter ATC Auto Temperature compensation Probe

With BNC cable connector and semipermeable membrane tip allows electrons to pass Dip in solution bottles

Cat NO	Type dia	Size Shape	Unit
AC-308-PT	Electrode combined +ATC pH meter		
AC-308-S47	Buffer solution capsule pH4 pH7 set		



Electrode Reference / Calomel for pH meter

With BNC cable connector and semipermeable membrane tip allows electrons to pass Dip in solution bottles

Cat NO	Type dia	Size Shape	Unit
AC-308-R	Electrode Reference for pH meter		
AC-308-PC	Electrode Calomel for pH meter		
AC-311-GI	Electrode Glass with banana Plug		



Electrode DP Probe Conductivity water test kit

With BNC cable connector and semipermeable membrane tip allows electrons to pass Dip in solution bottles

Cat NO	Type dia	Size Shape	Unit
AC-318-ATC	Electrode ATC Auto Temperature Probe		
AC-318-DO	Electrode DO Probe with membrane		
AC-318-CO	Electrode Conductivity Cell Probe Banana		



Electrode Glass reference Platinum for Potentiometer etc

With BNC cable / Banana Plug connector tip/ semipermeable membrane tip allows electrons to pass Dip in solution bottles

Cat NO	Type dia	Size Shape	Unit
AC-311-PT	Electrode Platinum with banana Plug		
AC-311-R	Electrode Reference with banana Plug		
AC-311-GI	Electrode Glass with banana Plug		
AC-311-SL	Electrode Silver with banana Plug		

