


Labomed Digital Camera iVu5100 5mp (Abron Cat no AB-37g51)



Dimensions (Units: mm)



| | |
|--|--|
| Sensor | 1/2.5" 5.0 MP Color CMOS |
| Active area | 5.70mm x 4.28mm |
| Sensor resolution | 2592 x 1944 pixels |
| Pixel Size | 2.2µm x 2.2µm |
| A/D conversion | 12 Bit RGB |
| Pixel clock | 48 MHz |
| Dynamic range | 68 dB (measured at 10 ms exposure) |
| Max. exposure | 3 seconds (Automatic Exp, WB, BL) |
| Frame rate (image size) | 30 fps (1080p) |
| Image resolution (standard binning) | 2592 x 1944 1280 x 720 |
| Interface | USB 2.0, HDMI, Secure Digital (SD) Card |
| Optical connection | C- Mount |
| Voltage supply | 5V/ 2 Amp (via USB 2.0 connector, optional external power adapter) |
| Power consumption | approx. 1.8W |
| Power indication | Green LED indicating power on and ready for capture |
| SD Card Slot | Accepts SD and SDHC (High Capacity) cards |
| Ambient conditions | Temperature: -30°C... +70°C Humidity: 5%... 80%, not condensing |
| Dimensions (LxWxH) | 145mm x 93mm x 123mm |
| Weight | 600g |
| Computer Requirements | PC: Microsoft Windows XP/ Vista/ 7 OS, 3 Ghz CPU, 2GB RAM, USB 2.0 |

| | |
|---|--|
| Package Contents <ul style="list-style-type: none"> • CapturePro Image Capture Software • Firewire IEEE 1394 Cable (3ft) | Warranty Information 2 Years Camera |
|---|--|