



Fluorometer Digital Photo 0-200 ppm AC-304F2

Highly sensitive, a full scale reading of 199.9
 Accuracy 0.2 ppm of quinine sulphate solution
 in 0.1 N Sulphuric acid
 Primary filter of peak transmission at 366 nm

Features:-

- 3-1/2 Digits Bright Red Seven Segment Led Display.
- High Stable & Accurate.
- Adjustable In 4 Ranges.0-0.2,2,2, 200(199) ppm
- Minimum Sample 5ml
- Determines The Concentration Of Fluorophores Quinine.
- Highly Sensitive Wide Range Photodiode

AC-304F2 Fluorimeter LED 0-200ppm With Mercury lamp primary filter peak at 366 nm



Abrón Exports
 www.abronexports.com

Fluorimeter digital LCD pen type with built in probe electrode AC-304FP

RANGE 0.1 to 10ppm x 01ppm
 ATC Temperature 0-60°Cx0.1°C
 Range of fluoride: 0 to 9.00 pF ±0.01 pF
 built-in automatic temperature compensation.
 easy to use for Field or laboratory testing.
 1 or 2 point calibration & adjustment.



AC-304FS Fluorimeter test strips with colour chart



AC-304H Fluorimeter Chemicals Test Kit with colour chart



Abrón Exports
 www.abronexports.com