

**Projection Abron Microscope Student
With screen 15cm dia Bulb AB-14**

Ideal for group teaching,
Focusing- for coarse/fine
Stage- 125x 125mm
Monocular Head- High Resolution
Lighting in base Halogen 6V 20w 220v AC
Mechanical Stage- 145x 125mm

Eye Piece	Achromatic Objective	Magnification
10X	10 X	100 X
	20 X	200 X

Cat No	Specifications	Type	Dia
AB-14-P	Projection Microscope with bulb screen Eye Piece 10x Objectives 10x 20x		15cm
AB-14-PR	Projection Microscope with bulb screen Eye Piece 10x Objectives 10x 20x+Mech stage		15cm



Projector Micro Slide halogen bulb +supply

Projection Distance 1 - 3.5 meter Resolution 1280 x 800
 Halogen lamp 24V 150 Watts Transformer 220V 50Hz AC

Cat No	Specifications	Type	Bulb
AB-45-P	Micro Slide Projector 15mm slide	3x1"	24v150W



**Senior Projection Microscope Big Screen
Halogen Bulb 12v 100w AB-46**

Body Steel sheet Powder coated
 screen 160mm / 200mm Rotatable
Lights Halogen Bulb 12V 100watt Halogen 220V AC
 Incident & Below sample light
 Objective Nose Piece Triple revolving nose piece
 Stage X-Y movement on ball bearing slides.
 Measuring By Micrometer (50mm 0.005mm
 Magnification Selectable



Abron Exports
www.abronexports.com

Cat No	Specifications	Magnification
AB-46-P6	Projection Microscope halogen screen Eye Piece 10X, Objective Magnification 50X, 100x 200X	160mm 5x 10x 20x
AB-46-P8	Projection Microscope halogen screen Objective 5x 10x 20x 40x Magnification 50X, 100x 200X ,400X With mechanical stage	200mm



AB-46-P8M	Projection Microscope halogen screen Objective 5x 10x 20x 40x,60x Magnification 100X, 200x 400X ,800X 1200x With mechanical stage	200mm
-----------	--	-------