

Ripple Tank Perspex Chamber AP-885PX

AP-885PX12

Ripple Tank Perspex tray chamber 12x12 inch
30cmx30cm



AP-885PXK Ripple Tank setup with supply
with Hand stroboscope



Ripple Tank Perspex Chamber Export Model AP-885P

Perspex frame water chamber with 4 legs, motor 1.5vdc with vibrating hanging by rubber for creating wave, bulb light for cell holder or power supply

a	Water chamber Perspex with 4 legs leveling	1 No
B	L Shape 75x25mm Barriers	4 no
C	L shape obstacles 15mm, 25mm 35mm wide 1 each	
D	C shape 6inch 15cm long circle	1 No
E	Mesh wave absorber triangular 20x25mm	4 No
F	Glass 100x50mm reflector	1 No
G	Motor with 2 balls for round waves	1 pair
H	Plan flat wave generator strip 100mm long	1 No
I	Paper white for putting on ground	1 No
J	Crocodile clips for holding Rubber bands	1 pair
K	Motor with Rubber bands set 1.5-3VDC	1 No
L	Motor Angle setting bridge	1 No
M	Bulb with holder shade wire plug 220v/12v	1 No
N	Stand arm for bulb holder	1 No
O	assembly and working Manual	1 No

OPTIONAL

AP-885P Ripple
Tank spare Power
supply for motor



AP-885S Ripple
Tank spare
Stroboscope hand
moving



AP-885SE Ripple
Tank spare Stroboscope
Electric



AP-885MP Ripple
Tank spare motor
vibrating strings



Abron Exports
www.abronexports.com