

Telephone No. 0091 171 2699478

Fax No: 0091 171 2698459

ABRON EXPORTS

New Add. :- 105 Industrial Estate HSIIDC, AMBALA CANTT. 133 001

E-mail No:- info@abronexports.com website <http://abronexports.com>

Probe Sonicator

Cat no : AC-386



Telephone No. 0091 171 2699478

Fax No: 0091 171 2698459

ABRON EXPORTS

New Add. :- 105 Industrial Estate HSIIDC, AMBALA CANTT. 133 001

E-mail No:- info@abronexports.com website <http://abronexports.com>



Probe Sonicator

SALIENT FEATURES

- Touch screen display - 4.3 inches TFT.
- User Password Protection available.
- Microprocessor controlled and completely programmable.
- Large LCD display with Time setting .
- Temperature Indicator and controller.
- Auto-tuning for convenience and optimal processing efficiency.
- Variable Power Output control.
- Auto frequency tracking.
- Output amplitude 0-100% adjustment.
- Automatic fault alarm.

SPECIFICATIONS

MODEL	AC-386
Frequency	20-25KHz - Automatic tracking
Process capacity	0.5-500ml (approx.)
Standard Probe (Titanium Alloy)	6mm (1/4") Standard
Display	4.3 inches TFT
Display Function	Temperature, Power, Time, etc.
Temperature range(°C)	0-99°C
Timer (Selectable)	Yes
Overload Protection	Yes
User Password Protection	Available
Input method & Data Memory	Touch screen control & 20 groups
Output Power	6.5-650W (adjustable)
Power supply	AC 220V /50HZ
Packing Dimension (LxWxH) & Gross Weight (Main Instrument)	480 x 325 x 450mm & 18 kg
Packing Dimension (LxWxH) & Gross Weight (Sound Proof Enclosure)	355 X 355 X 545mm & 10kg

* Specifications are subject to change without Notice

Standard Accessories: Main Instrument, Transducer with 6mm Probe, Temperature Probe, Power Cable, Soundproof Box Set, 3A Fuse, Instruction Manual, Test Certificate.

Optional Probe : 2mm, 3mm, 8mm, 10mm, 12mm